


DFN Distal Femoral Nail

Surgical Technique

Standard Locking

Spiral Blade Locking





511.790 Quick Coupling for Kirschner Wires

511.701 COMPACT™ AIR DRIVE II

511.760 Quick Coupling
for DHS/DCS® Triple Reamers

511.750 Quick Coupling



351.110 Reverse Awl



393.100 Universal Chuck with T-Handle




355.041 Guide Rod Ø 3.0 mm,
with flat tip

117.500 Reaming Instrument Set for
Medullary Nails (not illustrated)



399.430 Hammer 700 g



399.505 Hammer, synthetic



314.280 Holding Sleeve, large



357.360 Extraction Screw for
Spiral Blade



357.133 Extraction Screw for
AFN/DFN



357.117 Hammer Guide for DFN



357.026 Slotted Hammer 400 g,
can be mounted






Table of contents

Indications/contraindications	2
Implants	3
Standard locking	4
Spiral blade locking	17
Implant removal	20
Cleaning of instruments	23



Image intensifier control

Warning


This description is not sufficient for immediate application of the instrumentation. Instruction by a surgeon experienced in handling this instrumentation is highly recommended.

Indications/contraindications

The Distal Femoral Nail DFN is indicated for the stabilization of fractures of the distal femur. It can also be used for diaphyseal fractures in which a retrograde approach is indicated (e.g. ipsilateral tibia and/or patella fractures, proximal or distal endoprosthesis, adipositas permagna). These include according to the AO classification:

Indications

- Fractures of type 33-A1 to A3
- Fractures of type 33-C1 to C3.1
- Fractures of type 32-A to C



Contraindications

- Fractures of type 33-B, 33-C3.2 and 33-C3.3
- Proximal femoral fractures and high subtrochanteric fractures

Nails

Short and long, for standard locking and spiral blade locking
Solid: Ø 9.0 and 10.0 mm
Cannulated: Ø 12.0 mm

– Long Nail (451.705–711, 725–731, 825–831)

Lengths: 300, 340, 380, 420 mm

Proximal locking:

- 1 dynamic hole AP
- 1 static hole ML
- 1 static hole AP

– Short Nail (451.701–703, 721–723, 821–823)

Lengths: 160, 200, 240 mm

Proximal locking:

- 2 static holes ML

Proximal locking

Locking Bolts Ø 4.9 mm, self-tapping (459.260–680)

Lengths: 26 to 68 mm

Standard locking

For non-osteoporotic bone and simple supracondylar fracture morphology.

Locking Screws Ø 6.0 mm, self-tapping (450.862–872) (turquoise)

Lengths: 50 to 100 mm

End Cap for Locking Screw Ø 6.0 mm (451.896) (turquoise)

Spiral blade locking

For better purchase in osteoporotic bone and for all complex supracondylar fracture morphologies.

Locking Screws Ø 6.0 mm, self-tapping (450.862–872) (turquoise)

Lengths: 50 to 100 mm

Spiral Blades for DFN (450.882–892) (pink)

Lengths: 50 to 100 mm

End Cap for Spiral Blade (451.895) (pink)



Standard locking

Careful preoperative planning with clear classification of the fracture and correct choice of implants are mandatory for a successful surgical result.

1

Position patient

Position the patient supine on the operating table. The knee of the injured leg should be flexed 70°–90°. A leg roll may be used or the lower leg be flexed by lowering part of the table to allow proper reduction and stabilization of the reduced fracture. Position the image intensifier in such a way that visualisation of the entire femur is possible in anterior-posterior (AP) as well as lateral views.

2

Reduce fracture



If possible, carry out closed reduction under image intensifier control. In some cases the use of the Large Distractor (394.350) is required.

Note: After repositioning, and prior to the insertion of the nail, intra-articular fractures should be definitely stabilised with two interfragmentary lag screws. Position the screws so as not to interfere with the path of the nail. The intramedullary canal and the nail are slightly dorsal, therefore, lag screws are placed ventrally.


3

Determine nail length



Position the image intensifier for an AP view of the distal femur (Position 1). With long forceps, hold the Radiographic Ruler (357.122) parallel to and on the same level as the femur. Adjust the ruler until the distal end is at the desired nail insertion depth. Mark the skin laterally.

Move the image intensifier to the proximal femur (Position 2), place the distal end of the ruler on the skin mark, and take an AP image of the middle third of the femur (for short nails) and of the proximal femur (for long nails). Verify fracture reduction and read the nail length from the ruler image. When using a long nail, select its length so that it can be locked level with the lesser Trochanter.




4

Determine nail diameter



Place the radiographic ruler perpendicular to the femur and position it over the isthmus to confirm the nail diameter. Choose the nail diameter with which the transition between the medullary canal and the cortex is still visible left and right of the marking.



5

Incision

For 33-A.X and 32-X.X fractures, either make a transligamentar (ligamentum patellae) or a medial parapatellar incision. For 33-C.X fractures always make a medial parapatellar approach.



6

The Distal Femoral Nail may be used for either reamed or unreamed application.

Determine nail entry point and insert guide wire

The entry point for the nail is in the axis of the medullary canal and in the intercondylar notch, just anterior and lateral to the femoral attachment of the posterior cruciate ligament.

Insert the Protection Sleeve 17.0/15.0 (357.530) and the Drill Sleeve (357.531) through the incision to the bone. Insert the calibrated Guide Wire \varnothing 3.2 mm (357.129) to a depth of 10 to 15 cm, taking into consideration the 7° to 9° valgus angle of the anatomic axis of the femur. AP and lateral image intensifier control is mandatory.



7

Open medullary canal


Place the cannulated Drill Bit Ø 13.0 mm (351.270) into a cannulated power drill. Guide the drill bit through the protection sleeve and over the guide wire to the bone. Drill to a depth of approximately 30 mm. Remove the protection sleeve, the drill sleeve and the guide wire. Do not re-use the guide wire.

Remove the bone and cartilage debris and thoroughly rinse the knee joint.

Alternative

The medullary canal can also be opened with the Large Reverse Awl (351.110).

Reaming of the femur may be indicated depending on the fracture type or if the medullary canal is too narrow. If reaming is not indicated proceed with step 9.




8

Ream medullary cavity with SynReam (Option)

1 Assemble reduction system


Assemble the reduction system of the SynReam Intra-medullary Reaming System (189.060): Attach the T-Handle (351.150) at the rear of the SynReam Flexible Shaft (352.040) and a Reduction Head (352.050 or 352.055) at the front.



2 Reduce fracture


To secure the reduction head, insert the SynReam Reaming Rod \varnothing 2.5 mm (352.032, length 950 mm or 352.033, length 1150 mm) in a retrograde direction up to the olive. The olive must be located in the reduction head throughout reduction. Insert the assembled reduction system with the SynReam reaming rod into the medullary cavity and reduce the distal fragments under image intensifier control.

Note: Always reduce with the reaming rod, since secure fixation can only be ensured if the reduction system is used in conjunction with the reaming rod. Using the reduction system without the reaming rod entails the risk of losing the reduction head in the medullary canal.



3 Remove reduction instruments


After completing the reduction, remove the reduction instruments with the exception of the reaming rod, which must remain in the medullary cavity.



4 Assemble the reaming system


Connect the SynReam Flexible Shaft (352.040) to the drill and plug the first SynReam Medullary Reamer Head (352.085) onto the shaft. The reamer heads can be picked up directly, without hand contact, from the insert for medullary reamer heads using the SynReam flexible shaft.

Start with the smallest reamer head (\varnothing 8.5 mm, 352.085) and then increase in 0.5 mm increments using the larger reamer heads (352.090–190). The reaming depth should be identical to the chosen nail length.



5 Insert reaming system


Insert the assembled reaming system, without rotating it, over the SynReam reaming rod into the medullary canal. Use the Tissue Protector (351.050) to spare the soft tissues.



6 Ream medullary canal

Ream the medullary canal according to the standard procedure. Advance the reamer slowly and steadily at maximum drill speed. Secure the SynReam reaming rod with the Holding Forceps for SynReam Reaming Rod (351.782) to prevent it from rotating during reaming.


Note: Only ream over the SynReam Reaming Rod \varnothing 2.5 mm (352.032, length 950 mm or 352.033, length 1150 mm), since the rod ensures that a secure connection is maintained between the reamer and the flexible shaft.



7 Change reamer head

Having reamed the medullary cavity along its full length, retract the SynReam flexible shaft with the first reamer head until the whole reamer head is visible. Grasp the reaming rod with the Holding Forceps for SynReam Reaming Rod (351.782) immediately above the bone insertion point and hold in situ to avoid loss of reduction. Draw the SynReam flexible shaft through the slot of the Insert with Removing Device for SynReam Medullary Reamers (689.063) so as to remove the used reamer head without touching it.

The reamer head of the next size up can be picked up directly, without hand contact, from the insert for medullary reamer heads using the SynReam flexible shaft.



8 Complete medullary reaming

Repeat steps 5 to 7 for each additional reamer head (352.090–190) until the medullary canal is reamed to the desired diameter. Reaming is usually performed with increments of 0.5 mm.

Note: Only ream over the SynReam Reaming Rod \varnothing 2.5 mm (352.032 or 352.033), since the rod ensures that a secure connection is maintained between the reamer head and the flexible shaft.



Loosen blocked medullary reamer heads with left-right turns or with gentle hammer taps to the Holding Forceps for SynReam Reaming Rod (351.782) fastened to the SynReam Reaming Rod \varnothing 2.5 mm (352.032 or 352.033).

9

Attach nail to insertion handle

With the Insertion Handle (357.112) facing laterally, match the teeth on the handle with the notches in the nail. The anterior curvature of the implant must be aligned with that of the femur. Note the inscription ANTERIOR on the nail.

Slide the Connecting Screw (357.135) through the insertion handle and screw it into the nail. Tighten the connecting screw with the Pin Wrench (321.170) without stripping it.




10

Insert nail

If hammering with the Slotted Hammer (357.026) is necessary for insertion, attach the Hammer Guide (357.117) to the connecting screw and tighten the assembly.

Insert the nail by hand as far as possible twisting gently. When using a guide rod, slide the Cannulated Nail Ø 12.0 mm over the guide rod into the femoral canal. Monitor nail passage through the fracture under image intensification.


If necessary, insert the nail using light hammer blows with either the slotted hammer and hammer guide or with a Hammer (399.430/399.505). Using the distal end of the connecting screw as striking surface advance the nail until its distal end is 2 to 5 mm beyond the articular cartilage. The insertion depth is indicated by the grooves on the shaft of the insertion handle. Remove the guide rod.



Note: After nail insertion, check that the connecting screw is still tight as it could have loosened during hammering. Do not attach the aiming arm for standard locking or spiral blade locking, respectively, to the insertion handle until the nail is fully inserted. Otherwise, the aiming arm may loosen during nail insertion.

The standard locking is performed with two Locking Screws Ø 6.0 mm (450.862–872).


Spiral blade locking see page 17 (steps 11 to 17).



11

Mount aiming arm and insert trocar assembly


Attach the Aiming Arm for Standard Locking (357.115) to the insertion handle. Insert the turquoise triple trocar assembly (Protection Sleeve 11.0/8.0 (357.118), Drill Sleeve 8.0/4.9 (357.097) and Trocar Ø 4.9 mm (357.098)) through the cranial hole of the aiming arm and, after a stab incision, to the bone. Remove the trocar.



12


Drill and measure length for locking screw

Drill through both cortices using the calibrated Drill Bit Ø 4.9 mm (357.099) stopping just after penetrating the far cortex. Confirm the drill bit position radiographically. Make sure the drill sleeve is in contact with the bone and read the length for the locking screw from the ringmarks on the drill bit. The locking screw should not surpass the medial cortex by more than 2 mm.



Alternative


The Depth Gauge for Locking Bolts (357.791) may be used to determine the length of the locking screw. Remove the Drill Sleeve 8.0/4.9 and insert the measuring hook through the Protection Sleeve 11.0/8.0. Determine the length of the locking screw by adding 2 mm to the reading on the depth gauge.



13

Insert first locking screw

Insert the locking screw through the Protection Sleeve 11.0/8.0 using the Hexagonal Screwdriver (314.750) and make sure the screw head is in contact with the medial cortex. The locking screw should surpass the medial cortex by at least 2 mm, but not more than 4 mm.



14

Insert second locking screw

Repeat steps 11 to 13 for inserting the second locking screw in the caudal screw hole on the aiming arm.

Note: For standard locking use two locking screws whenever possible.



15

Proximal locking


Proximal locking for nails of 160–240 mm in length is done using two Locking Bolts Ø 4.9 mm which are inserted laterally. Nails of 300 mm or more are provided with a static hole in lateral direction, and a static hole and a dynamic locking slot in AP direction. These locking options make immediate or secondary dynamization of the fracture possible. The AP approach facilitates locking at the lesser trochanter level.

Align image


Check fracture reduction, the correct alignment of the fragments, and the correct length of the femur.



Align the image intensifier so that the distal hole is visible as a perfect circle in the centre of the image.



Round (correct)




Oval (incorrect)

16

Make incision



Determine the incision point and make a stab incision.



Determine
incision point

17

Drill




Under image intensification, insert the tip of the Drill Bit Ø 4.0 mm for the radiolucent drive through the incision and place it on the bone.


Tilt the drive until the drill bit appears centred in the locking hole. The drill bit will nearly fill out the locking hole image. Holding the drill in this position, drill through both cortices until the tip of the drill bit just penetrates the medial cortex.

Alternative

If there is no radiolucent drive available and locking is performed with the standard freehand technique use the three-flute Drill Bit Ø 4.0 mm (315.400).



Centre drill bit




Drill

18

Measure length for locking bolt and insert locking bolt

- ⓘ Measure the length for the locking bolt with the Depth Gauge for Locking Bolts (357.791). The length of the locking bolt is determined by adding at least 2, but not more than 4 mm to the length read from the depth gauge.
- ⓘ Insert the locking bolt using the hexagonal screwdriver.



19

Insert end cap


The use of the turquoise End Cap for Locking Screw (451.896) is mandatory. It fulfills two functions: It prevents bone ingrowth into the nail and blocks the distal Locking Screw \varnothing 6.0 mm thus providing a stable angle.

Remove the insertion instruments.

- ⓘ Using the hexagonal screwdriver, align the end cap with the nail axis and fix it by means of the holding sleeve. To minimise the chance of cross-threading, turn the end cap counterclockwise until the thread of the end cap aligns with that of the nail. By turning clockwise, completely screw the end cap into the nail, so that the locking screw is fixed. The last turns will offer some resistance caused by a groove in the thread which prevents the screw from loosening. Tighten the end cap firmly.

Alternative

The end cap can also be inserted after step 14 (page 14).



Spiral blade locking

The spiral blade locking option with one Spiral Blade (450.882–892) and one Locking Screw \varnothing 6.0 mm (450.862–872) is used for patients with an osteoporotic bone structure and/or extensive and complex supra-condylar comminuted fracture areas. The increased purchase area of the blade provides an optimized load distribution in the condyles and improves fixation of the distal fragment thus decreasing the risk of nail protrusion into the knee joint.


Perform steps 1 to 10 as described on pages 4 to 12.

11

Mount aiming arm and insert locking screw

Attach the Aiming Arm for Spiral Blade Locking (357.116) to the insertion handle. Insert the turquoise triple trocar assembly (Protection Sleeve 11.0/8.0 (357.118), Drill Sleeve 8.0/4.9 (357.097) and Trocar \varnothing 4.9 mm (357.098)) through the cranial hole of the aiming arm and, after a stab incision, to the bone. Remove the trocar.


Drilling, length measuring, and insertion of the Locking Screw \varnothing 6.0 mm are performed as described in steps 12 and 13 on pages 13 and 14.



12

Insert double-sleeve combination


Insert the pink double-sleeve combination (Protection Sleeve 15.0/13.0 (357.123) and Drill Sleeve 13.0/3.2 (357.124)) through the caudal hole on the aiming arm. Make a stab incision and advance the combination to the lateral cortex.



13

Insert guide wire and measure length for spiral blade


Insert the calibrated Guide Wire Ø 3.2 mm (357.129) into the condyles until its tip is flush with the medial cortex. Confirm wire position radiographically. Make sure the drill sleeve is in contact with the bone and read the length for the spiral blade from the guide wire. Remove the drill sleeve.



14

Open cortex


Insert the Drill Bit Ø 13.0 mm (351.270) over the guide wire through the protection sleeve and open the lateral cortex manually or with a power tool. The drill bit should open the lateral cortex only to prepare the seat of the spiral blade head. The automatic stop prevents the drill bit from penetrating too far. Remove the drill bit and protection sleeve.



15

Attach spiral blade to spiral inserter

Assemble the spiral blade insertion instruments using the Connecting Screw for Spiral Blades (357.132) and the Spiral Inserter (357.120). Attach the spiral blade of the previously determined length firmly to the spiral inserter.




16

Insert spiral blade

Bring the spiral blade and the spiral inserter over the guide wire through the aiming arm and to the bone. Insert the blade manually or by gentle hammer taps to the connecting screw attached to the spiral inserter. The correct insertion depth is reached when the blade head is flush with the lateral cortex. Verify under image intensification.

Detach the insertion instruments from the blade and remove the guide wire.



17

Proximal locking

Perform the proximal locking as described in steps 15 to 18 on pages 15 and 16.

18

Insert end cap


The use of the pink End Cap for Spiral Blade (451.895) is mandatory. It fulfills two functions: It prevents bone ingrowth into the nail and locks the spiral blade proximally providing a stable angle. The end cap presses against the edge of the spiral blade preventing it from backing out.

Remove the insertion instruments.

Using the hexagonal screwdriver, align the end cap with the nail axis and fix it by means of the holding sleeve. To minimise the chance of cross-threading, turn the end cap counterclockwise until the thread of the end cap aligns with that of the nail. By turning clockwise, completely screw the end cap into the nail, so that the spiral blade is fixed. The last turns will offer some resistance caused by a groove in the thread which prevents the screw from loosening. Tighten the end cap firmly.

Alternative

The end cap can also be inserted after step 16.




Implant removal

1

Remove end cap

Remove the ingrown tissue from the hexagonal recess of the end cap and the locking implants. Unscrew the end cap using the Hexagonal Screwdriver (314.750).




2

Remove the most distal locking implant

Remove the most distal locking implant (locking screw \varnothing 6 mm or spiral blade).

Locking screw \varnothing 6 mm: unscrew the locking screw using the Hexagonal Screwdriver (314.750) and Holding Sleeve (314.280).


Spiral blade: screw the Extraction Screw for Spiral Blade (357.360) into the blade head and screw the Hammer Guide (357.117) onto the extraction screw. Knock the spiral blade out with the Slide Hammer (357.026).



3

Remove locking bolts

Unscrew the locking bolts using the Hexagonal Screwdriver (314.750). If the bolts cannot be unscrewed because the thread no longer grips the bone, use the Holding Sleeve (314.280).



4


Mount hammer guide

Before removing the last locking implant, screw the Extraction Screw for DFN (357.133) into the distal end of the nail and tighten with the Pin Wrench (321.170). The last locking implant prevents the nail from rotating or sliding away.

Screw the Hammer Guide (357.117) onto the extraction screw for DFN and tighten.

Alternative


If the extraction screw for DFN is not available, the hammer guide can be screwed directly into the distal end of the nail.



5

Remove last locking implant


Unscrew the last locking implant using the Hexagonal Screwdriver (314.750) and, if necessary, the Holding Sleeve (314.280).



6

Remove nail

Knock the nail out with gentle blows of the slide hammer.



Cleaning of instruments

Cleaning of instruments

Proper functioning of the instruments can be maintained by careful cleaning.

DFN Distal Femoral Nail

Surgical Technique

Europe

Czech Republic

Stratec Medical spol. s r. o.
CZ-162 00 Praha 6
Tel. +420 2 2051 50 06
Fax +420 2 2051 50 06

Denmark

Stratec Medical A/S
DK-2730 Herlev
Tel. +45 44 53 45 44
Fax +45 44 53 27 09

Finland

Stratec Medical Oy
FIN-00180 Helsinki
Tel. +358 9 696 2540
Fax +358 9 6962 5420

Germany

Synthes GmbH
D-79224 Umkirch
Tel. +49 7665 50 30
Fax +49 7665 50 32 40

Great Britain

Stratec Medical Ltd
Welwyn Garden City
GB-Hertfordshire AL7 1LG
Tel. +44 1707 33 22 12
Fax +44 1707 33 85 04

Hungary

Stratec Medical Kft.
H-1021 Budapest
Tel. +36 1 391 45 60
Fax +36 1 391 45 69

Ireland

Fannin Healthcare
IRL-Dublin 18
Tel. +353 1 295 3401
Fax +353 1 295 4777

Norway

Stratec Medical A/S
N-0216 Oslo
Tel. +47 24 12 66 80
Fax +47 24 12 66 81

Poland

Stratec Medical Poland Sp. z o.o.
PL-02-222 Warszawa
Tel. +48 22 608 50 50
Fax +48 22 608 50 51

Portugal

Stratec Medical Lda.
P-1050-132 Lisboa
Tel. +351 21 357 00 54
Fax +351 21 357 68 39

Slovakia

Intes Poprad s.r.o.
SK-05899 Poprad
Tel. +421 52 7723 488
Fax +421 52 7723 393

Spain

Stratec Medical S.A.
E-28230 Las Rozas / Madrid
Tel. +34 902 190 128
Fax +34 902 190 347

Sweden

Stratec Medical AB
S-171 28 Solna
Tel. +46 8 743 71 00
Fax +46 8 641 48 56

Switzerland

Stratec Medical
CH-4436 Oberdorf
Tel. +41 61 965 61 11
Fax +41 61 965 66 00

Latin America

Argentina

Synthes Argentina S.A.
C1183AEY Buenos Aires
Tel. +54 11 4867 49 49
Fax +54 11 4867 49 55

Bolivia

Salur S.R.L.
La Paz
Tel. +591 2 223 390
Fax +591 2 223 390

Brazil

Synthes Ind. Com. Ltda.
CEP-13505-650 Rio Claro, SP
Tel. +55 19 3527 23 22
Fax +55 19 3527 04 82

Chile

Helico Artículos Médicos Ltda.
Santiago de Chile
Tel. +56 2 204 59 90
Fax +56 2 225 37 75

Costa Rica

Lorwen, S.A.
La Uruca, contiguo a Rapi Freno
San José
Tel. +506 290 5530
Fax +506 232 0511

Colombia

Synthes Colombia S.A.
Santafé de Bogotá D.C.
Tel. +57 1 612 11 20
Fax +57 1 612 979 35

Ecuador

Traumamed S.A.
Quito
Tel. +593 2 55 59 09
Fax +593 2 55 59 25

Mexico

Synthes S.M.P., S.A. de C.V.
53100 Ciudad Satélite
Tel. +52 5 562 40 59
Fax +52 5 562 11 85

Paraguay

Aldo M. Bergonzi Importaciones
Asunción
Tel. +595 21 29 23 00
Fax +595 21 29 23 00

Peru

Hoscli S.A.
Lima 34
Tel. +51 1 447 25 40
Fax +51 1 446 93 29

Uruguay

Resimpex S.A.
11200 Montevideo
Tel. +598 2 403 0677
Fax +598 2 401 2868

Venezuela

I.P.M. Ingeniería y Productos Médicos C.A.
Caracas
Tel. +58 2 577 10 44
Fax +58 2 574 22 57