

# Blood Gas RAPIDPoint 500e

## Operator Competency Assessment/Access Application

Name (please print):		Employee ID #:	
Department (state your Department regardless of whether there is a blood gas machine there. This is for management purposes):		Hospital Email:	
Blood Gas Access barcode required? (i.e Are you a new user, or is your barcode damaged/lost?)		<input type="checkbox"/> Yes	<input type="checkbox"/> No
Training Type	<input type="checkbox"/> Initial/New	<input type="checkbox"/> Renewal of Certification	<input type="checkbox"/> Renewal of Competency
(please tick relevant competency)			
Blood Gas Machine Site Access Required:	<input type="checkbox"/> RCH Emergency	<input type="checkbox"/> RCH PICU	<input type="checkbox"/> RCH NICU
	<input type="checkbox"/> RWH NICU	<input type="checkbox"/> RWH 5 South	<input type="checkbox"/> RCH Theatres

For certification on this device, the Operator must demonstrate the following competencies:	Yes	N/A
<b>Pre-analytcs</b>		
Identifies components of the RapidPoint 500e System		
States the patient identification requirements for this test type		
Describes the sample types and requirements (e.g arterial/venous/capillary/mixed venous/panels {eg. GLU only, Hb only}, well mixed, blood heparin ratio, minimum sample volume.)		
Can describe the effect of the following concerns on the sample: <ul style="list-style-type: none"> <li>- Clotted specimen (blotting sample on gauze pre-testing)</li> <li>- Diluted specimen</li> <li>- Settled red cells (mixing technique)</li> <li>- Delayed analysis (stable at Room temperature for 10 minutes when kept mixed)</li> <li>- Haemolysis</li> <li>- Cleaning product contamination (sodium interference)</li> </ul>		
Describe the process for removal of air bubbles prior to mixing and before initiating testing		
<b>Analytcs</b>		
States that Operator Identification must NOT be shared and able to log onto the RapidPoint 500e analyser		
Recognise and interpret the analyser status (i.e calibrating, sensor offline)		
Selects the correct sample type		
Ensures Positive Patient Identification prior to beginning the test		
Correctly introduces sample to the analyser		
<b>Post-Analytcs</b>		
Removes and correctly disposes of sample (i.e yellow clinical waste)		
Understands the printed results and escalates critical results to the appropriate treating team		
Can access and reprint Patient Results		
<b>Maintenance</b>		
Can describe clot management <ul style="list-style-type: none"> <li>- Obstructions/Insufficient/Bubbles in sample</li> <li>- Describe how to replace the sample port</li> </ul>		
Describes procedure for cleaning analyser with 0.5% bleach (supplied by the laboratory) and bench area/barcode reader		
Can describe where to find device procedure		
<b>Operator Statement:</b>		
By signing this form, I declare that I have completed the training required to a competent level to run blood gas samples and have been approved by the Nurse Unit Manager/Trainer to gain a personal login for the RapidPoint 500e. I understand that on-going competency will be required at defined intervals.		
Operator Signature:	Date:	
Trainer's Name (please print):	Trainer's Signature:	
<input type="checkbox"/> Unit Staff	<input type="checkbox"/> Laboratory Staff	<input type="checkbox"/> Vendor
	Date:	
<b>LAB USE ONLY</b>	Profile <input type="checkbox"/>	Barcode <input type="checkbox"/>
	Certified	Yrs
Comments:		