

Cobe Spectra™ General Maintenance Record

Refer to: Chapter 10 General Maintenance Cobe Spectra™ Apheresis System. Essentials Guide. Gambro BCT. Lakewood Colorado, AP-E-001, for a detailed description of completion of the set tasks. Please note, weekly and monthly maintenance is completed at the discretion of the operator and in accordance with the number of procedures performed (eg, maintenance not required in weeks machine is not used). Please sign as verification of completion of the set task.

MONTH:

| TASK | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | 25 | 26 | 27 | 28 | 29 | 30 | 31 |
|----------------------------|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| POST PROCEDURE | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| WEEKLY | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Inlet Air Detector | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Return Air Detector | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| RBC Detector | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Fluid Leak Detector | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Centrifuge Pressure Sensor | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Clean Door Tracks | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| MONTHLY | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Front Door Sensor | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Pressure Pod Seals | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| CARIDIANBCT | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |