

Thermometer Check Record

WORKING THERMOMETER (WT) IDENTIFICATION NUMBER: _____

| DATE | RT ID | VISUAL INSPECTION ACCEPTABLE (Y/N) | BATTERY CHANGED (Y/N/NA) | TEMP. CHECKED (°C) | WT TEMP. (°C) | RT TEMP. (°C) | CORRECTED RT TEMP. (°C) <i>(CHECK Test Certificate)</i> | WT-RT _(corrected) (°C) | CORRECTION REQUIRED (°C) | SIGNED |
|------|-------|------------------------------------|--------------------------|--------------------|---------------|---------------|--|-----------------------------------|--------------------------|--------|
| | | RT: | RT: | IP | | | | | | |
| | | | | 4°C | | | | | | |
| | | WT: | WT: | AT | | | | | | |
| | | RT: | RT: | IP | | | | | | |
| | | | | 4°C | | | | | | |
| | | WT: | WT: | AT | | | | | | |
| | | RT: | RT: | IP | | | | | | |
| | | | | 4°C | | | | | | |
| | | WT: | WT: | AT | | | | | | |
| | | RT: | RT: | IP | | | | | | |
| | | | | 4°C | | | | | | |
| | | WT: | WT: | AT | | | | | | |

(REFERENCE THERMOMETER (RT) managed by Laboratory Services, IP=ice-point, AT=ambient temperature)

NOTE: WT should be within ±1°C of the RT. If there is more than 1°C, consult the Quality Manager.

SUMMARY

| DATE | ACTION IF APPLICABLE | AUTHORISATION |
|------|----------------------|---------------|
| | | |
| | | |
| | | |

REFERENCE DOCUMENT: Thermometer Checks. QS-W-001.