











Nasal Cannula







F&P Optiflow Junior				
Product	Item Code	Approx. Weight Range	Max. Flow Rate (L/Min)	Spare Wigglepads
Optiflow Junior Nasal Cannula				
 Premature Size	OPT312	<2kg	8	OPT010
 Neonatal Size	OPT314	1-8kg	8	OPT012
 Infant Size	OPT316	3-15kg	20	OPT012
 Pediatric Size	OPT318	12-22kg	25	OPT012
Optiflow Junior Tubing Kit	RT330 for MR850, 900PT531 for AIRVO 2	N/A	Refer to Optiflow Junior Nasal	N/A

1 APPLY CANNULA



2 REMOVE CANNULA



3 REPLACE WIGGLEPADS™

