

COVID-19 PIPER airway equipment



Covidien
viral filter
<10kg



Pharma
mini
viral filter
>10kg



BMV/T piece/Neopuff



Tapes/
NGT/
Hollister



Oropharyngeal
airways

**Pre/Re/
Apnoeic
oxygenation**

Video
laryngoscopy
with correct size
blade

Age	ETT size (Microcuff)	ETT size (Mallinckrodt)	Blade size	Depth (Oral)
2.5 – 3.5 kg	3.0 uncuffed		0 – 1	8.5 – 9.5 cm
3.5 kg – 8 mo	3.0 cuffed		0 – 1	9.5 – 10.5 cm
8 mo – 2 yrs	3.5		1	11 cm
2 – 4 yrs	4		1 – 2	12 cm
4 – 6 yrs	4.5		2	14 cm
6 – 8 yrs	5		2	15 cm
8 – 10 yrs	5.5		2 – 3	16 cm
10 – 12 yrs	6		2 – 3	17 cm
12 – 14 yrs		6.5	3	18 cm
>14 yrs		7	3	19 cm

Magill
forceps


Clamps

Stylet

Laryngoscope #2

Bougie

ETT size below



Emesis bag
for laryngoscope

In-line
suction



Plan B

LMA

Weight (kg)	LMA size
Neonates: IGEL	
1.5 – 5	1
5 – 12	2
Paediatrics: Ambu	
<5kg	1
5 – 10k	1.5
10 – 20	2
20 – 30	2
30 – 50	2.5
50 – 70	3
70 – 100	4
12 – 14 yrs	5
>14 yrs	